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Reagent kit for immuno-turbidimetric determination of Complement C4 in human serum.

Principle

Measurement of antigen/antibody reaction by the endpoint method.

Reference values

Serum: 9- 36 mg/dl

It is recommended that each laboratory should assign its own normal range.

Reagents

- 1. Reagent 1 (R1)**
buffer
sodium-azide (0,50 g/l)
- 2. Reagent 2 (R2)**
antiserum
sodium-azide (0,50 g/l)

Precaution

Discard cloudy reagent. Avoid contamination by using clean laboratory material (pipettes, plastic vials for analyzers,...). The reagents contain sodium azide. To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent.

Stability of the reagents

Store the reagent protected from light. Store it between 2-8 °C.
 without opening: till the expiry date indicated on the label
 after opening: 14 days
 calibration frequency: 7 days
 onboard stability: 7-14 days
 Stability data are valid only when using new system bottle!

Sample

Use fresh serum. The serum can be stored at 2-8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

PROCEDURE

The reagents are ready for use.

Assay conditions

Wavelength: 340 nm
 Temperature: 37 °C
 Cuvette: 1 cm light path
 Read against: sample blank
 Method: endpoint (increasing)

Pipette into cuvette

	Blank	Standard	Sample
1. reagent (R1)	360 µl	360 µl	360 µl
Standard	3 µl or 3 µl	3 µl	
Sample			3 µl

Mix, wait 1 minute then add:

2. reagent (R2)		90 µl	90 µl
------------------------	--	-------	-------

Mix and after 7 minutes incubation read the absorbance against the sample blank.

Calibration:

5 level (at the least 3) standard set is necessary for calibration.

S1: Distilled water or physiological salt solution
 S2-S6: 5 level standard set

New calibration is recommended:
 - after reagent lot change,
 - as required following quality control procedures.

Calculation

Based on the calibration curve.

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits.

PERFORMANCES DATA

The following data were obtained using the Olympus AU600 and Advia 1650 analysers (37°C).

Linearity

The test is linear up to 120 mg/dl (1,2 g/l).

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. The estimated value of limit of detection (LOD) is 0 mg/dl.

Precision

	Reproducibility		
	Average concentration (g/l)	SD	CV %
sample I	0,33	0,03	9,30
sample II	0,31	0,02	0,76
sample III	0,34	0,03	8,74
sample IV	0,29	0,02	5,33

Sample	Repeatability		
	Average concentration (g/l)	SD	CV%
sample I.	0,12	0,004	3,02
sample II.	0,27	0,006	2,05
sample III.	0,32	0,003	0,99

Correlation

Comparative studies were done to compare our reagent with another commercial Complement C4 reagent. The results from these studies are detailed below.
 Correlation coefficient: $r=0.9719$
 Linear regression: $y (g/l)= 0,910x + 0,009$
 (x= other commercial reagent, y= own reagent).

Specificity

Bilirubin 20 mg/dl, triglycerides 2500 mg/dl, hemoglobin 1000 mg/dl don't interfere with the assay up to the given levels.

NOTE

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

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