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Verzió: 42-DL-2015-04

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**Reagent kit for immuno-turbidimetric determination of Immunoglobulin M (IgM) in human serum.**

**Principle**

Measurement of antigen/antibody reaction by the endpoint method.

**Reference values**

**Serum: 70-250 mg/dl**

It is recommended that each laboratory should assign its own normal range.

**Reagents**

**1. Reagent 1 (R1)**

buffer  
sodium-azide (0,50 g/l)

**2. Reagent 2 (R2)**

antiserum  
sodium-azide (0,50g/l)

**Precaution**

Discard cloudy reagent. Avoid contamination by using clean laboratory material (pipettes, plastic vials for analyzers,...). The reagents contain sodium azide. To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent.

**Stability of the reagents**

Store the reagents protected from light. Store it between 2-8 °C.  
 without opening: till the expiry date indicated on the label  
 after opening: 14 days  
 calibration frequency: 7 days  
 onboard stability: 7-14 days  
 Stability data are valid only when using new system bottle!

**Sample**

Use fresh serum. The serum can be stored at 2-8 °C for 48 hours. If stored for a longer period, the sample should be frozen.

**PROCEDURE**

The reagents are ready for use.

**Assay conditions**

Wavelength: 340 nm  
 Temperature: 37 °C  
 Cuvette: 1 cm light path  
 Read against: sample blank  
 Method: endpoint (increasing)

**Pipette into cuvette**

	Blank	Standard	Sample
<b>1. reagent (R1)</b>	320 µl	320 µl	320 µl
<b>Standard</b>	3 µl or	3 µl	
<b>Sample</b>	3 µl		3 µl

Mix, wait 20 seconds then add:

<b>2. reagent (R2)</b>		80 µl	80 µl
------------------------	--	-------	-------

Mix and after 4 minutes incubation read the absorbance against the sample blank.

**Calibration:**

5 level (at the least 3) standard set for calibration is necessary.

S1: Distilled water or physiological salt solution

S2-S6: 5 level standard set

New calibration is recommended:

- after reagent lot change,
- as required following quality control procedures.

**Calculation**

Based on the calibration curve.

**Quality control**

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

**PERFORMANCES DATA**

The following data were obtained using the Olympus AU600 and Advia 1650 analysers (37°C).

**Linearity**

The test is linear up to 500 mg/dl ( 5 g/l).

**Sensitivity**

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. The estimated value of limit of detection (LOD) is 0 mg/dl.

**Precision**

	Reproducibility		
	Average concentration (g/l)	SD	CV %
sample I.	1,63	0,16	9,71
sample II.	1,56	0,12	7,92
sample III.	1,39	0,15	10,52
sample IV	1,27	0,08	6,26

	Repeatability		
	Average concentration (g/l)	SD	CV%
sample I.	1,59	0,01	0,71
sample II.	1,15	0,01	0,77
sample III.	8,34	0,09	1,06

**Correlation**

Comparative studies were done to compare our reagent with another commercial IgM reagent. The results from these studies are detailed below.  
 Correlation coefficient:  $r=0,9981$   
 Linear regression:  $y (g/l)= 0,966x + 0,078$   
 (x= other commercial reagent, y= own reagent).

**Specificity**


Bilirubin 20 mg/dl, triglycerides 2500 mg/dl, haemoglobin 1000 mg/dl don't interfere with the assay up to the given levels.

**NOTE**

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

**For in vitro diagnostic use only.**

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

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