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Reagent kit for determination of albumin concentration in serum. Colorimetric bromocresol green method.

Albumin is a protein that is formed in the liver. Approximately 50-60 % of total protein is albumin, the rest is globulin. The albumin helps to maintain normal distribution of water in the body (colloidal osmotic pressure), and also helps in the transport of blood constituents such as ions, pigments, bilirubin, hormones, enzymes, and drugs. The ratio of albumin to globulin is an important indicator of certain disease states. Albumin level is decreased in several liver diseases, malabsorption, diarrhoea, eclampsia and nephrosis. The level is increased in dehydration.

Principle

At pH=4.2, albumin bind with bromocresol green to produce a blue-green complex. The change in absorbance at 628 nm correlates with the concentration of albumin.

Reference values

Serum albumin 35-50 g/l (3,5-5 g/dl)

It is recommended that each laboratory should assign its own normal range.

Reagents

1. Reagent (R1)

bromocresol green	0,15 mmol/l
succinate buffer pH=4.20	75 mmol/l
Brij 35	7 ml/l

Samples

Serum free of haemolysis. Stability: 20 days (2-8°C)

Stability

Avoid light and store at 2-8°C

without opening:	till the date indicated on the label
after opening:	30 days
calibration frequency:	7 days
onboard stability:	7-30 days

Stability data are valid only when using new system bottle!

PROCEDURE

Working reagent

The reagent is ready for use.

If the absorbance of working reagent is higher than 0.2 at 578 nm the reagent can not be used.

Assay conditions

Wavelength:	628 (578-650) nm
Temperature:	37 °C
Cuvette:	1 cm light path
Read against:	reagent blank
Method:	endpoint (increasing)

Pipette into cuvette

	Blank	Standard	
Reagent	1 ml	1 ml	1 ml
Distilled water	10 µl		
Standard		10 µl	
Sample			10 µl

Mix and read the optical density (A) after a 3-minute incubation.

Calibration: (37°C, Bromocresol green method)

S1: Distilled water
S2: Diagnosticum DunaCal Cat. No.: Deal or Albumin standard Cat. No.: 123456 or Roche C.F.A.S. (Calibrator for automated system)
Randox Calibration Serum Level I or Randox Calibration Serum Level II

Calibration frequency

Two-point calibration is recommended:
- after reagent lot change,
- as required following quality control procedures.

Calculation

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance
C = Concentration

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish

corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

PERFORMANCE DATA

The following data were obtained using the Olympus 600 analyser (37°C). Conversion factor: [g/l]=10×[g/dl]

Linearity

The test is linear between 10 - 69 g/l (1 - 6,9 g/dl)

Limit of detection

The limit of detection is 0,187 g/l (0,0187 g/dl)

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 0.06 g/l (0,006g/dl) albumin concentration at 578 nm.

Precision

n=20	Average concentration (g/l)	SD	CV%
sample I	34.2	0.41	1.19
sample II	32.3	0.52	1.60

Correlation

Comparative studies were done to compare our reagent with another commercial albumin reagent. The results from these studies are detailed below.

Correlation coefficient: r=0.9969

Linear regression: y (g/l)= 0.976x+0.653

(x= other commercial reagent, y= own reagent).

Specificity


Bilirubin 855 µmol/l (50mg/dl), lipid 1000 mg/dl, glucose 55.5 mmol/l (1000mg/dl) and ascorbic acid 2.84 mmol/l (50mg/dl) don't interfere with the assay up to the given levels.

NOTE

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

Bibliography

Doumas B. et al.: Clin. chem. Acta 31, 87 (1971)
Doumas B. et al.: In Standard Methods of Clinical Chemistry, Acad. Press N.Y., 7, 175 (1972)
Drupt F.: Pharm. Biol. 9, 777 (1977)
Szabó A.: Klinikai laboratóriumi vizsgálatok és paraméterek (2010) (ISBN 978-963-9879-75-1)

Version: 12-DL-2016-12
Date of revision: 2016-12