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Verzió: 16-DL-2016-13

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

Reagent kit for determination of chloride ion concentration in serum, urine and cerebrospinal fluid. A colorimetric endpoint method based on the reaction with mercuric thiocyanate.

Chloride ion represents that anion of the salt and water household which is present in the highest concentrations in the body. It is primarily found as NaCl in the extracellular compartment and as HCl in the gastrointestinal tract. Low chloride concentrations are associated with severe vomiting, diarrhoea, colitis ulcerosa, diabetic acidosis, Addison's disease. Decreased chloride concentrations are found also in cases when drugs that need chloride ions for their absorption are taken for prolonged time. Increased levels are observed in cases of dehydration, congestive heart failure, Cushing's syndrome, hyperventilation, anaemia, nephritis and renal obstruction.

Principle

Chloride ion in acidic environment in presence of ferric nitrate forms a colored complex with mercuric thiocyanate. Intensity of the developed colour is proportional to the chloride ion concentration in the sample.

Reference values

Serum: 98-107 mmol/l (348-380 mg/dl)
Cerebrospinal fluid: 119-131 mmol/l (421-465 mg/dl)
Urine: 110-250 mmol/24h (3,89-8,87 g/24h)
27-371 mmol/l (95,9-1317 mg/dl)

It is recommended that each laboratory should assign its own normal range.

Reagents

1. Reagent (R1)

Mercuric thiocyanate 2 mmol/l
Ferric nitrate 20 mmol/l
Nitric acid 29 mmol/l

Samples

Serum, stability: 1 week (2-8°). Urine (1+2 dilution), cerebrospinal fluid.

Stability

without opening: till the expiry date indicated on the label
after opening: 28 days
calibration frequency: 7 days
onboard stability: 7-28 days
Stability data are valid only when using new system bottle!

PROCEDURE

Reagent is ready for use.

If the absorbance of working reagent is higher than 0.1 at 492 nm the reagent can not be used.

Assay conditions

Wavelength: 500 (480-520) nm
Temperature: 37°C
Cuvette: 1 cm light path
Method: endpoint (increasing)

Pipette into cuvette

| | Blank | Standard | Sample |
|-----------------|-------|----------|--------|
| Reagent | 1 ml | 1 ml | 1 ml |
| Distilled water | 10 µl | | |
| Standard | | 10 µl | |
| Sample | | | 10 µl |

Mix and read the absorbance against the blank.

Calibration: (37°C, mercuric thiocyanate method)

S1: Distilled water
S2: Diagnosticum DunaCal Cat. No.: Deacal
Chloride standard Cat. no.: 351301 or
Roche C.F.A.S. (Calibrator for automated system)
Randox Calibration Serum Level I or
Randox Calibration Serum Level II

Calibration frequency

Two point calibration is recommended:
- after reagent lot change,
- as required following quality control procedures.

Calculation using calibration

$$\frac{A_{\text{sample}}}{A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = absorbance, C = concentration

Quality control

Quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P). DunaCont Urine two levels control (DCONU2) is recommended for urine analysis.

PERFORMANCES DATA

The following data were obtained using the Hitachi 717 analyzer (37°C). Conversion factor: [mmol/l]=[mg/dl]/3,55

Linearity

Between 14-135 mmol/l (49,7-479 mg/dl)

Limit of detection

The limit of detection is 0,27 mmol/l (0,96 mg/dl)

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 0.25 mmol/l (0,89 mg/dl) chloride concentration at 492 nm.

Precision

Reproducibility

| n=20 | serum | | urine | | | |
|------|---------------|------|-------|---------------|------|-----|
| | sample mmol/l | SD | CV% | sample mmol/l | SD | CV% |
| 111 | 2,6 | 2,35 | 59,56 | 1,39 | 2,33 | |
| 109 | 1,3 | 1,19 | 172 | 4,6 | 2,67 | |

Repeatability

| n=20 | serum | | |
|------|--------------|------|-----|
| | serum mmol/l | SD | CV% |
| 83,4 | 1,19 | 1,42 | |
| 112 | 1,05 | 0,94 | |

Correlation

Comparative studies were done to compare our reagent with ion-selective measurement. The results from these studies are detailed below.

Correlation coefficient: r=0.9903

Linear regression: y (µmol/l)= 0.999x-4.366

(x= ion selective measurement, y= own reagent).

Specificity

Hemoglobin 77.5 µmol/l (500mg/dl), bilirubin 855 µmol/l (50 mg/dl), lipid 700mg/dl, glucose 55.5 mmol/l (1000mg/dl) and ascorbic acid 1.99 mmol/l (35mg/dl) don't interfere with the assay up to the given levels.

NOTE

In the course of determination please use disposable equipments. Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

| | | | |
|---|----------------------------|---|------------------|
|  | In vitro diagnostic device |  | Batch code |
|  | Manufacturer |  | Catalogue number |
|  | CE-marking |  | This way up |
|  | Temperature limitations |  | Biological risk |
|  | Use by (year/month) | | |

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