| | | | |
|-------------|---------------------------------|------------------------------------|---------------------------------|
| Kódszám: | 41751,41751S | 41752,41752S | 41753 |
| Kiszerezés: | 1x500 ml (1x250 ml+1x250 ml) | 10x500 ml (10x250 ml+10x250 ml) | 10x40 ml (10x20 ml+10x20 ml) |

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Reagent kit for the determination of creatinine concentration in serum and urine. A colorimetric, alkaline picrate method (Jaffé).

Creatinine is released during metabolism of creatine phosphate, and is excreted by the kidneys. Creatinine concentration in blood and in urine represents a primary indicator for renal function, especially that for glomerular filtration. Increased levels are associated with acute renal impairment, chronic nephritis, obstruction of the urinary tract, strong physical overloading. Low creatinine concentrations are found in conditions with juvenile diabetes mellitus, pregnancy and muscular dystrophy.

Principle

Creatinine forms with alkaline picrate (in ratio of 1:1) a colored creatinine picrate complex containing ionic bonds. The rate of formation of the colored complex is proportional to the creatinine concentration.

Reference values

| | | |
|---------------|----------------|---|
| Serum: | Male: | 60-110 µmol/l (0,60-1,2 mg/dl) |
| | Female: | 50-100 µmol/l (0,50-1,10 mg/dl) |
| Urine: | Male: | 2,12-3,46 mmol/l (24-392 mg/dl) |
| | Female: | 1,41-2,89 mmol/l (16-327 mg/dl) |
| | | 5-18 mmol/24 hours (600 - 2000 mg/24h) |

It is recommended that each laboratory should assign its own normal range.

Reagents

1. Reagent (R1)

| | |
|--------------------|------------|
| Sodium hydroxyde | 400 mmol/l |
| Disodium phosphate | 25 mmol/l |

2. Reagent (R2)

| | |
|-------------|-------------|
| Picric acid | 17.5 mmol/l |
|-------------|-------------|

Safety instructions:

Reagent 1:

Warning. Contains Sodium-hydroxide. Causes skin irritation. Causes serious eye irritation. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if present and easy to do – continue rinsing.

Reagent 2:

EUH001: Explosive when dry

Samples

Serum and 12 h or 24 h collected urine, resp. Urine must be diluted in ratio of 1:49 with distilled water. Stability in serum: 1 week (2-8°C)

Stability of the reagents

| | |
|---|---|
| without opening: | till the expiry date indicated on the label |
| after opening: | 30 days |
| calibration frequency: | 1 days |
| onboard stability: | 1-3 days |
| Stability data are valid only when using new system bottle! | |

PROCEDURE

Preparation of working reagent

Mix (R1) with (R2) in a ratio of 1:1.

Stability of working reagent

at: 20-25°C 1 month

If the absorbance of working reagent is higher than 0.4 at 492 nm the reagent can not be used.

Assay conditions

| | |
|--------------|----------------------|
| Wavelength: | 492 (480-520) nm |
| Temperature: | 37 °C |
| Cuvette: | 1 cm light path |
| Method: | kinetic (increasing) |

Pipette into cuvette

| | Standard | Sample |
|------------------------|----------|--------|
| Standard | 100 µl | |
| Sample | | 100 µl |
| Working reagent | 1 ml | 1 ml |

Mix and after 30 seconds read the absorbance against distilled water (A1). After 2 minutes incubation read the absorbance again (A2).

The reagent kit is suitable for two-reagent method, too. Reagents (R1) and (R2) can also be pipetted separately (0.5-0.5 ml).

Calibration: (37°C, Jaffé method)

- S1: Distilled water
- S2: Diagnosticum DunaCal Cat. No.: Dcal or Creatinine standard Cat. No.: 950911 or Roche C.F.A.S. (Calibrator for automated system)
- Randox Calibration Serum Level II

Calibration frequency

- Two point calibration is recommended:
- after reagent lot change,
- as required following quality control procedures.

Calculation using calibration

$$\frac{A2_{sample} - A1_{sample}}{A2_{standard} - A1_{standard}} \times C_{standard} = C_{sample}$$

A = absorbance
 C = concentration

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P). DunaCont Urine two levels control (DCONU2) is recommended for urine analysis.

PERFORMANCES DATA

The following data were obtained using the Olympus 600 analyzer (37°C). Conversion factor: [µmol/l]=[mg/dl]×88,42

Linearity

Between 20 - 1326 µmol/l (0,23 - 15 mg/dl).

Limit of detection: 1,89 µmol/l (0,02 mg/dl)

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs units/min is equivalent to 3.82 µmol/l (0,04 mg/dl) creatinine concentration at 492 nm.

Precision

Reproducibility

| n=20 | serum | | urine | | |
|------|---------------|------|-------|---------------|------|
| | sample µmol/l | SD | CV% | sample µmol/l | SD |
| 110 | 0,98 | 0,88 | 6406 | 586 | 9,16 |
| 338 | 3,36 | 0,99 | 7130 | 750 | 10,5 |

Repeatability

| n=20 | serum | |
|------|---------------|------|
| | sample µmol/l | SD |
| 78.4 | 0,75 | 0,96 |
| 236 | 8,98 | 3,81 |

Correlation

Comparative studies were done to compare our reagent with another commercial Creatinine reagent.

The results from these studies are detailed below.

Correlation coefficient: r=0.9994 Linear regression: y (µmol/l)= 0.987x+8.523 (x= other commercial reagent, y= own reagent).

Specificity

Hemoglobin 6,2 µmol/l (40 mg/dl), bilirubin 68,4 µmol/l (4 mg/dl), lipid 1000mg/dl, glucose 2,8 mmol/l (50mg/dl) and ascorbic acid 0,85 mmol/l (15mg/dl) don't interfere with the assay up to the given levels.

Note

The presence of high concentrations of keton derivatives, pyruvic acid, ascorbic acid, glucose, urea, proteins interfere with the test. Kinetic assay reduces the disturbing effects of some substances. Because of differences in reaction rates, it is advisable to take readings during the first few minutes of the reaction. Creatinine reacts faster with picrate as does with some non-specific disturbing chromogens. Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

| | | | |
|---|----------------------------|---|------------------|
|  | In vitro diagnostic device |  | Batch code |
|  | Manufacturer |  | Catalogue number |
|  | CE-marking |  | This way up |
|  | Temperature limitations |  | Biological risk |
|  | Use by (year/month) | | |

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