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N. Watanabe N. et. a.: Clin. Chem. 32: 1551. (1986)

Verzió: 32-DL-2016-13

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Reagent kit for the quantitative determination of total protein concentration in urine and liquor. Pyrogallol Red, direct colorimetric method.

Principle

Protein molecules in urine or liquor bind to the molybdate pyrogallol Red complex. The formation of the protein-dye complex causes a shift in the absorbance maximum from 467 nm to 598 nm.

Reference values

Urine total protein <140 mg/24h
 Male: 10-160 mg/l (1-16 mg/dl)
 Female: 0-160 mg/l (0-16 mg/dl)
 Cerebrospinal liquor 80-320 mg/l (8,0-32,0 mg/dl)
 It is recommended that each laboratory should assign its own normal range.

Reagents

1. Color reagent (R1)

Pyrogallol-red 60 µmol/l
 Sodium molybdate 40 µmol/l
 Sodium oxalate 1.04 mmol/l
 Sodium benzoate 3.47 mmol/l
 Methanol 100 mol/l
 Succinic acid 50 mmol/l

Safety instructions:

Reagent 1:

Danger. Contains Methanol. May cause damage to organs. Do not breathe dust/fumes/gas/mist/vapours/spray.

IF SWALLOWED: Immediately call a POISON CENTER/doctor

Sample

Urine. Cerebrospinal fluid.

Stability of the reagent

Store the reagent protected from light. Store it between 2-8 °C.
 without opening: till the expiry date indicated on the label
 after opening: 14 days
 calibration frequency: 5 days
 onboard stability: 5-14 days
 Stability data are valid only when using new system bottle!

PROCEDURE

Working reagent

The reagents are ready for use. Keep away from light.
 If the absorbance of working reagent is higher than 0.25 at 578 nm the reagent can not be used.

Assay conditions

Wavelength: 598 nm
 Temperature: 37 °C
 Cuvette: 1 cm light path
 Read against: reagent blank
 Method: endpoint (increasing)

Pipette into cuvette

	Blank	Standard	Sample
Color reagent	1 ml	1 ml	1 ml
Sodium chloride (9 g/l)	20µl		
Standard		20µl	
Sample			20µl

Mix and incubate for 10 minutes at room temperature.

Calibration: (37°C, Pyrogallol Red method)

Two-point:
 S1: Distilled water
 S2: Total protein ultrasensitive standard (50,0 mg/dl) Cat. No.: 520715H
 Four-point:
 S1: Distilled water
 S2: Total protein ultrasensitive standard (200,0 mg/dl) Cat. No.: 520715N
 ratio 1:3 dilution (50,0 mg/dl)
 S3: Total protein ultrasensitive standard
 ratio 1:1 dilution (100,0 mg/dl)
 S4: Total protein ultrasensitive standard
 without dilution (200,0 mg/dl)
 Use distilled water for dilution!

Calibration frequency

Two point calibration is recommended:
 - after reagent lot change,
 - as required following quality control procedures.

Calculation using calibration

$$\frac{A_{sample}}{A_{standard}} \times C_{standard} = C_{sample}$$

A = Absorbance
 C = Concentration

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. DunaCont Urine two levels control (DCONU2) is recommended for urine analysis.

PERFORMANCES DATA

The following data were obtained using the Olympus 600 analyzer (37°C). Conversion factor: [mg/l]=[mg/dl]×10

Linearity

The test is linear up to 2000 mg/l (200 mg/dl) if four-point calibration is used, and 1500 mg/l (150 mg/dl) if the calibration is two-point. Lower limit of linearity is 50 mg/l (5 mg/dl)

Limit of detection 5 mg/l

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 1.9 mg/l (0,19 mg/dl) protein concentration at 578 nm.

Selectivity

Bilirubin up to 1 mg/dl (5 mg/dl protein) and 5 mg/dl (20 mg/dl protein) don't interfere with the assay up to the given levels.

Precision

Reproducibility

n=20	serum		urine		
sample mg/l	SD	CV%	sample mg/l	SD	CV%
114	2,26	1,98	75	2,3	3,11
727	20,5	2,82	537	18,5	3,45

Repeatability

n=20	serum	
sample mg/l	SD	CV%
26,6	1,23	4,61
2448	12,9	0,53

Correlation

Comparative studies were done to compare our reagent with another commercial Total protein reagent.

The results from these studies are detailed below.

Correlation coefficient: $r = 0,9778$

Linear regression: $y \text{ (mg/l)} = 0,993x - 1,657$

(x= other commercial reagent, y= own reagent).

NOTE

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

Bibliography

N. Watanabe N. et. a.: Clin. Chem. 32: 1551. (1986)

Version: 32-DL-2016-13

Date of revision: 2016-12