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Verzió: 36-DL-2017-11

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Reagent kit for the quantitative determination of direct bilirubin in serum. Diazo-sulfanilic acid method.

Principle

Sulfanilic acid reacts with sodium nitrite to form diazotized sulfanilic acid. In the presence of dimethylsulfoxide, total bilirubin reacts with diazotized sulfanilic acid to form azobilirubin. In the absence of dimethylsulfoxide, only the direct bilirubin reacts to give azobilirubin.

Reference value

Serum: direct bilirubin <5 µmol/l (<0,3mg/dl)
 total bilirubin 3-17 µmol/l (0,2-1 mg/dl)

It is recommended that each laboratory should assign its own normal range.

Reagents

Direct Bilirubin

1. Reagent (RD1)
 Sulfanilic acid 3,2 mmol/l
 Hydrochloric acid 165 mmol/l

2. Reagent (R2)
 Sodium nitrite 8.6 mmol/l

Sample

Serum free of haemolysis. Stability 7 days at 2-8°C.

Bilirubin in serum is light sensitive and it is recommended that serum be stored in the dark. The bilirubin decomposes fast in the serum and controls.

The reagent is not suitable for bilirubin determination of infants.

Procedure

All reagents are ready for use.

Avoid direct exposure to light!

Preparation of working reagent

Mixing ratio of RD1 and R2: 125 ml/25 ml

Reagents can only be applied with previous mixing!

Stability

without opening: till the expiry date indicated on the label
 after opening: 30 days
 calibration frequency: 1 days
 onboard stability: 1-3 days
 Stability data are valid only when using new system bottle!

Stability of the working reagent

If the absorbancy of working reagent is higher than 0,02 at 546 nm the reagent can not be used.

Assay Conditions

Wavelength: 555 nm (540-560 nm)
 Secondary wavelength: 600 nm
 Temperature: 37 °C
 Cuvette: 1 cm light path
 Measure: end point
 Read against: reagent blank

Pipette into cuvette

	Sample	Calibrator
Working reagent	1,0 ml	1,0 ml
Sample	75 µl	-
Calibrator	-	75 µl

Mix and read the absorbancy (A) after exactly 3 minutes incubation.

Calibration: (37°C, Diazo-Sulfanilic Acid)

S1: Distilled water
 S2: Diagnosticum DunaCal Cat. No.: Dcal or Randox Calibration Serum Level I or II or Roche C.F.A.S. (Calibrator for automated system)

Calibration is recommended:

- after reagent lot change,
- as required following quality control procedures.

Calculation using calibration

$$\frac{A_{Sample}}{A_{Calibrator}} \times C_{Calibrator} = C_{Sample}$$

A = Absorbance, C = Concentration

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

PERFORMANCES DATA

The following data were obtained using the Olympus 600 analyzer. Conversion factor: [µmol/l]=17,1×[mg/dl]

Linearity

The test is linear between 1 - 100 µmol/l (0,06 - 5,88 mg/dl) bilirubin concentration (37°C).

Limit of detection 0,07 µmol/l (0,004 mg/dl)

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 0.575 µmol/l (0,03 mg/dl) Bilirubin concentration at 540 nm.

Precision

n=20	Reproducibility		
	Average concentration (µmol/l)	SD	CV%
sample I	13.8	0.282	2.04
sample II	37.8	1.053	2.79

Correlation

Comparative studies were done to compare our reagent with another commercial Bilirubin Direct reagent.

The results from these studies are detailed below.

Correlation coefficient: r=0.9957

Linear regression: y (µmol/l)= 1.05x-0.52

(x= other commercial reagent, y= own reagent).

Specificity


Lipid 400 mg/dl, glucose 55.5 mmol/l (1000 mg/dl) and ascorbic acid 2.84 mmol/l (50 mg/dl) don't interfere with the assay at the given levels.

Note

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein

For in vitro diagnostic use only!

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

Bibliography

Hijmans Van den Bergh A. A., Muller P.: *Biochem*, 77; 90, (1916).
 Walters M. I., Gerarde R. W.: *Microchem.*, 15; 231 (1970).
 Szabó A.: *Klinikai laboratóriumi vizsgálatok és paraméterek (2010) (ISBN 978-963-9879-75-1)*

Version: 36-DL-2017-11

Date of revision: 2017-01