
Kódszám:	46163	46161,46161S	46162,46162S
Kiszerezés:	10x20 ml (10x15 ml+ 10x5 ml)	120 ml (1x90 ml+ 1x30 ml)	500 ml (1x375 ml+1x125 ml)

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
### Reagent kit for determination of creatinine concentration in serum and urine. Colorimetric, enzymatic test.

Creatinine is released during metabolism of creatine phosphate, and is excreted by the kidneys. Creatinine concentration in blood and in urine represents a primary indicator for renal function, especially that for glomerular filtration. Increased levels are associated with acute renal impairment, chronic nephritis, obstruction of the urinary tract, strong physical overloading. Low creatinine concentrations are found in conditions with juvenile diabetes mellitus, pregnancy and muscular dystrophy.

#### Principle

#### Reference values

<b>Serum:</b>	<b>Male:</b>	<b>60-110 µmol/l (0,60-1,2 mg/dl)</b>
	<b>Female:</b>	<b>50-100 µmol/l (0,50-1,10 mg/dl)</b>
<b>Urine:</b>	<b>Male:</b>	<b>2,12-3,46 mmol/l (24-392 mg/dl)</b>
	<b>Female:</b>	<b>1,41-2,89 mmol/l (16-327 mg/dl)</b>
		<b>5-18 mmol/24 hours (600 - 2000 mg/24h)</b>

It is recommended that each laboratory should assign its own normal range.

#### Reagents

##### 1. Reagent (R1)

Creatinase	>20 kU/l
Sarcosine oxidase	>6 kU/l
Ascorbate oxidase	>2 kU/l
Catalase	>100 kU/l
TOOS	>0,40 mmol/l

##### 2. Reagent (R2)

4-aminoantipyrine	>2,5 mmol/l
Creatininase	>250 kU/l
Peroxidase	>50 kU/l

#### Precaution

Discard cloudy reagent. Avoid contamination by using clean laboratory materials (pipettes, plastic vials, ...) for analyzers.

#### Samples

Serum free of haemolysis. Stability: 1 week (2-25°C). Urine diluted in ratio of 1:49 with distilled water.

#### Stability of the reagents

Store the reagent protected from light. Store it between 2-8 °C.  
 without opening: till the expiry date indicated on the label  
 after opening: 28 days  
 calibration frequency: 7 days  
 onboard stability: 7-28 days  
 Stability data are valid only when using new system bottle!

### PROCEDURE

#### Working reagent

The reagents are ready for use. If the absorbance of reagent 1 is higher than 0.01 or the absorbance of the reagent 2 is higher than 0.015 at 546 nm the reagent can not be used.

#### Assay conditions

Wavelength:	555 (540-570) nm
Temperature:	37°C
Cuvette:	1 cm light path
Read against:	against reagent blank
Method:	endpoint (increasing)

#### Pipette into cuvette

	Blank	Standard	Sample
<b>R1</b>	1 ml	1 ml	1 ml
<b>Distilled water</b>	45µl		
<b>Standard</b>		45µl	
<b>Sample</b>			45µl

Mix and read the absorbance (A1) after a 5-minute incubation then add:

<b>R2</b>	330 µl	330 µl	330 µl
-----------	--------	--------	--------

Mix and read the absorbance (A2) after a 5-minute incubation.

#### Calibration: (37°C, enzymatic, colorimetric method)

S1: Distilled water  
 S2: Diagnosticum DunaCal Cat. No.: Dcal or Creatinine standard Cat. No.: 950911 or Roche C.F.A.S. (Calibrator for automated system)  
 Randox Calibration Serum Level I or

Randox Calibration Serum Level II

#### Calibration frequency:

Two point calibration is recommended  
 - after reagent lot change,  
 - as required following quality control procedures.

#### Calculation

$$\frac{(A2 - A1)_{\text{sample}}}{(A2 - A1)_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance  
 C = Concentration

#### Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P). DunaCont Urine two levels control (DCONU2) is recommended for urine analysis.

### PERFORMANCE DATA

The following data were obtained using the Olympus 600 analyzer (37°C). Conversion factor: [µmol/l]=[mg/dl]×88,42

#### Linearity

The test is linear up to 1770 µmol/l (20 mg/dl).

#### Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 2.2 µmol/l (0,025mg/dl) creatinine concentration at 546 nm.

#### Precision

n=20	Reproducibility		
	Average concentration (µmol/l)	SD	CV%
<b>Sample I</b>	112	0.94	0.84
<b>Sample II.</b>	352	3.48	0.99

n=20	Repeatability		
	Average concentration (µmol/l)	SD	CV%
<b>Sample I</b>	99	1.08	1.09
<b>Sample II.</b>	157	1.31	0.83

#### Correlation

Comparative studies were done to compare our reagent with another commercial Creatinine Jaffe reagent.

The results from these studies are detailed below.

Correlation coefficient: r = 0.9988  
 Linear regression: y (µmol/l) = 1.05x - 14.60  
 (x = other commercial reagent, y = own reagent).

#### Specificity

Bilirubin 684 µmol/l (40 mg/dl), ascorbic acid 200 mg/dl (11,36 mmol/l) and haemoglobin 80 µmol/l (500 mg/dl) don't interfere with the assay up to the given levels.

#### Note

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

#### For in vitro diagnostic use only.

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

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