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
### Reagent kit for the determination of the aspartate aminotransferase (AST) activity in serum based upon IFCC recommendations.

AST originates in various tissues and is a dimer molecule containing one molecule of Pyridoxal phosphate (coenzyme) in each monomer, which is essential to its catalytic activity. Depending on the sites of origin inside the cell there are two isoenzymes with different pH optimum: the mitochondrial m-AST, and the soluble cytosolic s-AST. The two isoenzymes can be separated by electrophoresis. The enzyme catalyses the transfer of amino groups during the metabolism of Amino acids and, alpha-Ketoacids. The activity of AST in the serum is significantly increased during heart, liver, kidney and muscle diseases (tissue injuries, functional disorders). The activity of the enzyme is increased 4-8 hours following a myocardial infarction, reaching its peak in 2-3 days and declining on the fifth and sixth days.

#### Principle

Two substrates participate in the reaction catalyzed by AST, L-aspartate and Oxaloglutarate. With the help of NADH coenzyme, Malate dehydrogenase (MDH) contained in the reagent catalyses the transformation of Oxalacetate released in the first reaction. The oxido-reductive process of NADH/NAD<sup>+</sup> is indicated by a decrease in absorbance at 340 nm. The Lactate dehydrogenase (LDH) in the medium counteracts the disturbing effect of Pyruvate contained in the sample.

#### Reference values

|         |                                    |
|---------|------------------------------------|
| Male:   | 5-45 U/l (0,085-0,765 $\mu$ kat/l) |
| Female: | 5-40 U/l (0,085-0,67 $\mu$ kat/l)  |

It is recommended that each laboratory should assign its own normal range.

#### Reagents

##### 1.Reagent (R1)

|                      |             |
|----------------------|-------------|
| Tris buffer, pH:7.80 | 88 mmol/l   |
| L-Aspartate          | 260 mmol/l  |
| LDH                  | 1500 U/l    |
| MDH                  | 900 U/l     |
| NADH                 | 0,24 mmol/l |

##### 2.Reagent (R2)

|                         |           |
|-------------------------|-----------|
| $\alpha$ -Ketoglutarate | 12 mmol/l |
|-------------------------|-----------|

#### Samples

Serum free of haemolysis. Haemolysis, lipaemia interferes with the test. Stability: 7 days (2-8°C)

#### Stability

|                        |   |
|------------------------|---|
| without opening:       | till the expiry date indicated on the label |
| after opening:         | 28 days                                     |
| calibration frequency: | 7 days                                      |
| onboard stability:     | 7-28 days                                   |

Stability data are valid only when using new system bottle!

#### Precaution

These reagents contain 0.1 % sodium azide. Discard cloudy reagent. To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent. Avoid contamination by using clean laboratory materials (pipettes, plastic vials for analyzers, ...).

### PROCEDURE

#### Preparation and stability of working reagent

- One-reagent procedure:

Mix 2 volumes of reagent 1 with 1 volume of reagent 2.

Stability: Keep the reagent protected from light.

|               |         |
|---------------|---------|
| at 20-25 °C : | 5 days  |
| at 2-8 °C:    | 4 weeks |

- Two-reagent procedure: reagents are ready for use.

If the absorbance of working reagent is lower than 1.2 at 334 nm the reagent can not be used.

#### Assay conditions

|               |                      |
|---------------|----------------------|
| Wavelength :  | 340 (334-365) nm     |
| Temperature : | 37°C                 |
| Cuvette :     | 1 cm light path      |
| Read against: | distilled water      |
| Method:       | kinetic (decreasing) |

#### One reagent procedure

|                   |             |
|-------------------|-------------|
| Working reagent   | 1 ml        |
| Sample or Control | 100 $\mu$ l |

Mix and after a 1-minute incubation, measure the change of absorbance per minute ( $\Delta A/\text{min}$ ) during 3 minutes.

#### Two reagent procedure

|                   |             |
|-------------------|-------------|
| R1                | 1 ml        |
| Sample or Control | 150 $\mu$ l |

Mix, incubate for one minute 37°C and add:

|    |             |
|----|-------------|
| R2 | 500 $\mu$ l |
|----|-------------|

Mix and after a 1-minute incubation, measure the change of absorbance per minute ( $\Delta A/\text{min}$ ) during 3 minutes.

#### Calibration (37°C, IFCC method without pyridoxal-phosphate)

S1: Distilled water

S2: Diagnosticum DunaCal Cat. No.: Deal or Roche C.F.A.S. (Calibrator for automated system) or Randox Calibration Serum Level I or Randox Calibration Serum Level II

#### Calibration frequency

Two point calibration is recommended

- after reagent lot change,

- as required following quality control procedures.

#### Calculation using calibration

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance

C = Concentration

#### Calculation using factor

340 nm: Activity (U/l)=  $\Delta A/\text{min}$ . x 2000; 334 nm: Activity (U/l)=  $\Delta A/\text{min}$ . x 1790

340 nm: Activity ( $\mu$ kat/l)=  $\Delta A/\text{min}$ . x 33,33; 334 nm: Activity ( $\mu$ kat/l)=  $\Delta A/\text{min}$ . x 54,24

#### Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

### PERFORMANCE DATA

The following data were obtained using the Olympus 400 analyzer (37°C). Conversion factor: [U/l]=[ $\mu$ kat/l]>60

#### Linearity

The method is linear in the range 5 – 260 U/l (0,083 – 4,33  $\mu$ kat/l)

#### Limit of detection

The limit of detection is 0,67 U/l

#### Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs units/min is equivalent to 1.79 U/l (0,03 $\mu$ kat/l) GOT activity at 334 nm.

#### Precision

| n=20<br>Activity U/l | repeatability |      | reproducibility |      |      |
|----------------------|---------------|------|-----------------|------|------|
|                      | SD            | CV%  | Activity U/l    | SD   | CV%  |
| 38,3                 | 0,79          | 2,07 | 19,1            | 0,19 | 1,02 |
| 134                  | 1,49          | 1,11 | 128             | 1,41 | 1,10 |
|                      |               |      | 227             | 1,55 | 0,68 |

#### Correlation

Comparative studies were done to compare our reagent with another commercial GOT reagent. The results from these studies are detailed below.

Correlation coefficient: r = 0.9998

Linear regression: y (U/l)= 0.992x+0.729

(x= other commercial reagent, y= own reagent).

#### Specificity

Hemoglobin 1.6  $\mu$ mol/l (10 mg/dl), bilirubin 855  $\mu$ mol/l (50 mg/dl), lipid 300 mg/dl, glucose 55.5 mmol/l (1000 mg/dl) and ascorbic acid 2.84 mmol/l (50 mg/dl) don't interfere with the assay up to the given levels.


#### Note

The test doesn't contain pyridoxal-phosphate.

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

#### For in vitro diagnostic use only.

#### The following symbols can be used on the labels

|   |                            |   |                  |
|---|----------------------------|---|------------------|
|  | In vitro diagnostic device |  | Batch code       |
|  | Manufacturer               |  | Catalogue number |
|  | CE-marking                 |  | This way up      |
|  | Temperature limitations    |  | Biological risk  |
|  | Use by (year/month)        |   |                  |

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