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

**Reagent kit for the determination of the alanine aminotransferase (ALT) activity in serum based upon IFCC recommendations.**

ALT is found in the cytosol of cells, it is a non-tissue specific soluble enzyme. ALT catalyses the transfer of amino groups during the transformations of aminoacids and alpha-ketoacids. Pyridoxal phosphate activates the process. The enzyme found in the serum is principally derived from the liver and kidney. The serum enzyme activity is increased during various hepatic disease states including hepatitis.

**Principle**

ALT catalyses the transformation of L-Alanin and 2-Oxoglutarate at optimal pH. The Pyruvate released in the reaction is transformed by Lactate dehydrogenase (LDH) in the presence of NADH /NAD<sup>+</sup> coenzyme to L-lactate, while the NADH/NAD<sup>+</sup> oxidoreductive process shows a decrease in absorbance at 340 nm. The change in absorbance correlates with serum ALT activity.

**Reference values**

Male	5-45 U/l (0,085-0,765 µkat/l)
Female	5-40 U/l (0,085-0,67 µkat/l)

It is recommended that each laboratory should assign its own normal range.

**Reagents**

**1.Reagent (R1)**

Tris buffer, pH:7.50	110 mmol/l
L-Alanine	600 mmol/l
LDH	1500 U/l
NADH	240 µmol/l

**2.Reagent (R2)**

α-Ketoglutarate	16 mmol/l
-----------------	-----------

**Precaution**

Discard cloudy reagent. These reagents contain 0.1% sodium azide. Avoid contamination by using clean laboratory materials (pipettes, plastic vials for analyzers, ...). To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent.

**Samples**

Serum free of haemolysis. Haemolysis, lipaemia interfere with the test. Stability: 1 week (2-8°C)

**Stability**

without opening:	till the expiry date indicated on the label
after opening:	30 days
calibration frequency:	3 days
onboard stability:	3-14 days

Stability data are valid only when using new system bottle!

**PROCEDURE**

**Preparation and stability of working reagent**

• One-reagent procedure:

Mix 2 volumes of reagent 1 with 1 volume of reagent 2.  
 Stability: at 20-25 °C: 5 days  
 at 2-8 °C: 2 weeks

• Two-reagent procedure: reagents are ready for use.

If the absorbance of working reagent is lower than 1.1 at 334 nm the reagent can not be used.

**Assay conditions**

Wavelength :	340 (334-365) nm
Temperature :	37°C
Cuvette :	1 cm light path
Read against:	distilled water
Method:	kinetic (decreasing)

• **One-reagent procedure**

Working reagent	1 ml
Sample or Control	100 µl

Mix and after a 2-minute incubation, measure the change of absorbance per minute (ΔA/min) during 2 minutes.

• **Two-reagent procedure**

R1	1 ml
Sample or Control	150 µl

Mix, incubate for one minute 37 °C and add:

R2	500 µl
----	--------

Mix and after a 2-minute incubation, measure the change of absorbance per minute (ΔA/min) during 2 minutes.

**Calibration**(37°C, IFCC without pyridoxal-phosphate)

S1: Distilled water  
 S2: Diagnosticum DunaCal Cat. No.: Deal or Roche C.F.A.S. (Calibrator for automated system) or

Randox Calibration Serum Level I or  
 Randox Calibration Serum Level II

**Calibration frequency**

Two point calibration is recommended  
 - after reagent lot change,  
 - as required following quality control procedures.

**Calculation using calibration**

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance, C = Concentration

**Calculation using factor**

340 nm: Activity (U/l) = ΔA/min. x 2200; (µkat/l) = ΔA/min. x 36,67  
 334 nm: Activity (U/l) = ΔA/min. x 1950; (µkat/l) = ΔA/min. x 32,50

**Quality control**

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

**PERFORMANCES DATA**

The following data were obtained using the Olympus 400 analyzer (37°C). Conversion factor: [U/l]=[µkat/l]×60

**Linearity**

The method is linear in the range 5 – 450 U/l (0,08 – 7,5 µkat/l)

**Limit of detection**

The limit of detection is 2,0 U/l

**Sensitivity**

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs units/min is equivalent to 1.95 U/l (0,033 µkat/l) GPT activity at 334 nm.

**Precision**

Reproducibility, n=20	Average activity (U/l)	SD	CV%
sample I	34.2	0.80	2.35
sample II	136	0.80	0.59

Repeatability, n=20	Average activity (U/l)	SD	CV%
sample I	14.5	0.158	1.09
sample II	318	3.032	0.95

**Correlation**

Comparative studies were done to compare our reagent with another commercial GPT reagent.

The results from these studies are detailed below.

Correlation coefficient: r=0.9997  
 Linear regression: y (U/l)= 1.001x+0.199  
 (x= other commercial reagent, y= own reagent).

**Specificity**

Hemoglobin 1.6 µmol/l (10 mg/dl), bilirubin 257 µmol/l (15 mg/dl), lipid 300 mg/dl, glucose 55.5 mmol/l (1000 mg/dl) and ascorbic acid 2.84 mmol/l (50mg/dl) don't interfere with the assay up to the given levels.

**Note**

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

**For in vitro diagnostic use only.**

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

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