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
**Kinetic determination of the lactate dehydrogenase activity in serum based upon DGKC recommendations.**

Lactate dehydrogenase (LDH) is present in every cell, it is a tetramere molecule which is a combination of two different tissue components (M-muscle, H-heart). There are five different isoenzymes: LDH-1: LDH-2: LDH-3: LDH-4: LDH-5 = 20: 34: 23: 12: 11. The serum activity is mainly composed of LDH-1, LDH-2 derived from the myocardium and red blood cells, and LDH-5 derived from the liver. The activities of isoenzymes are different in cases of certain substrates. The inhibitors and pH sensitivities are different. The various fractions were determined using chromatography in the past but more recently electrophoresis is the method of choice. The ratio of isoenzymes indicates certain disease states. The enzyme activity significantly increases 8-12 hours following a myocardial infarction and declines after 4-5 days. There is an increase in liver diseases, in certain anaemia and tissue injuries. The enzyme catalyses the Pyruvate / Lactate transformation at optimal pH.

**Principle**

LDH catalyses the transformation of Pyruvate to Lactate in pH=7.5 Tris buffer with NaCl in the presence of NADH coenzyme. The transformation of NADH to NAD<sup>+</sup> is accompanied by a decrease in absorbance at 340 nm. The change in absorbance correlates with the LDH activity in the serum.

**Reference values**

**Serum LDH-P activity: 240-480 U/l (4,0-8,0 µkat/l)**

It is recommended that each laboratory should assign its own normal range.

**Reagents**
**1.Reagent (R1)**

Phosphate buffer, pH=7.50 56 mmol/l  
 Pyruvate 1.6 mmol/l

**2.Reagent (R2)**

NADH 240 µmol/l

**Precautions**

Discard cloudy reagent. Avoid contamination by using clean laboratory material (pipettes, plastic vials for analyzers,...). These reagents contain sodium azide (0.1 %). To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent.

**Sample**

Serum free of haemolysis. Stability: 7 days (2-8°C)

**Stability**

without opening: till the expiry date indicated on the label  
 after opening: 35 days  
 calibration frequency: 7 days  
 onboard stability: 7-35 days  
 Stability data are valid only when using new system bottle!

**PROCEDURE**
**Preparation and stability of working reagent**
**• One-reagent procedure:**

Mix 1 volume of reagent 1 (R1) with 1 volume of reagent 2 (R2).

Stability : at 20-25°C: 1 week  
 at 2-8 °C: 4 weeks

**• Two-reagent procedure:**

The reagents are ready for use.

If the absorbance of working reagent is lower than 1.1 at 334 nm the reagent can not be used.

**Assay conditions**

Wavelength: 340 (334-365) nm  
 Temperature: 37 °C  
 Cuvette: 1 cm light path  
 Read against: distilled water  
 Method: kinetic (decreasing)

**• One-reagent procedure:**

Working reagent	1 ml
Sample	25µl

Mix and after 1 minute incubation, measure the change of optical density per minute (ΔA/min) during 3 minutes.

**• Two-reagent procedure:**

R1	1 ml
Sample	50µl

Mix and after 1 minute incubation, add:

R2	1 ml
----	------

Mix and after 1 minute incubation, measure the change of absorbance per minute (ΔA/min) during 3 minutes.

**Calibration:** (37°C, DGKC method)

S1: Distilled water  
 S2: Diagnosticum DunaCal Cat. No.: Dcal or Roche C.F.A.S. (Calibrator for automated system) or Randox Calibration Serum Level I  
 Calibration frequency  
 Two point calibration is recommended:

- after reagent lot change,  
 - as required following quality control procedures.

**Calculation using calibration**

$$\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance,  
 C = Concentration

**Calculation using factor**

334 nm: Activity (U/l) = ΔA/min x 8252; (µkat/l) = ΔA/min x 137  
 340 nm: Activity (U/l) = ΔA/min x 8450; (µkat/l) = ΔA/min x 141

**Quality control**

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

**PERFORMANCES DATA**

The following data were obtained using the Olympus 600 analyzer (37°C). Conversion factor: [U/l]=[µkat/l]x60

**Linearity**

The method is linear between 10 - 1200 U/l (20 µkat/l)

**Limit of detection 3,5 U/l (0,05 µkat/l)**
**Sensitivity**

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs units/min is equivalent to 8.252 U/l (0,138µkat/l) LDH-P activity at 334 nm.

**Precision**

n=20	Reproduceability		Repeatability		
sample U/l	SD	CV%	sample U/l	SD	CV%
344,6	3,41	0,99	244	21,5	8,9
523,6	17,13	3,27	435	24,2	5,6

**Correlation**

Comparative studies were done to compare our reagent with another commercial LDH-P reagent.

The results from these studies are detailed below.

Correlation coefficient: r = 0,9987

Linear regression: y (U/l) = 1.042x+13.1

(x= other commercial reagent, y= own reagent).

**Specificity**

Bilirubin 855µmol/l (50mg/dl), lipid 1000mg/dl, glucose 55.5 mmol/l (1000mg/dl) and ascorbic acid 2.84 mmol/l (50mg/dl) don't interfere with the assay up to the given levels.

**Note**

Haemolysis interferes with the test.

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

**For in vitro diagnostic use only.**
**The following symbols can be used on the labels**

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

**Bibliography**

Ann. Biol. Clin. 1982; 40:123.

Version: 64-DL-2017-01

Date of revision: 2017-01