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
Reagent kit for the quantitative determination of glucose concentration in serum, liquor and urine. Enzymatic test.

Determination of glucose concentration is important in the diagnosis and treatment of disorders of carbohydrate metabolism. Values higher or lower than the reference are of diagnostic significance. The levels are increased in diabetes mellitus, hyperthyroidism and in the hyperactivity of the pituitary gland. Decreased levels are observed in cases of overproduction of insulin by the pancreas, with tumors of the pancreas, as well as with hypofunction of the organs involved in glucose synthesis and carbohydrate metabolism.

Principle

HK = Hexokinase

G-6-P = Glucose-6-phosphate

G-6-PDH = Glucose-6-phosphate dehydrogenase

Reference values

Serum, adult:	3,6-6 mmol/l (64,5 – 107,8 mg/dl)
Serum, child:	3,3-5,6 mmol/l (60-100 mg/dl)
Cerebrospinal fluid:	2,78-3,89 mmol/l (50-69,9 mg/dl)
Urine:	0-1.1 mmol/l (0-19,8 mg/dl) <2,78 mmol/day (<500 mg/day)

It is recommended that each laboratory should assign its own normal range.

Reagents

1. Reagent (R1)

Pipes buffer, pH=7.60	80 mmol/l
NAD	3 mmol/l
ATP	1.7 mmol/l

2. Reagent (2)

Magnesium salt	4 mmol/l
Hexokinase	1700 U/l
G-6-PDH	1700 U/l

Precaution

Discard cloudy reagent. Avoid contamination by using clean laboratory material (pipettes, plastic vials for analyzers, ...). These reagents contain 0.1 % sodium azide. To avoid the possible build-up of azide compounds, flush waste-pipes with water after the disposal of undiluted reagent.

Sample

Serum free of haemolysis, urine, cerebrospinal fluid. Stability in serum: 72 hours (2-8°C)

Stability

without opening:	till the expiry date indicated on the label
after opening:	30 days
calibration frequency:	7 days
onboard stability:	7-30 days

Stability data are valid only when using new system bottle!

PROCEDURE

Preparation and stability of working reagent

- One-reagent procedure

Mix 4 volumes of the reagent 1 (R1) with 1 volume of the reagent 2 (R2).

Stability:	at 20-25°C: 2 weeks
	at 2-8°C: 1 month

Store the reagent protected from light.

- Two-reagent procedure

Reagents are ready to use.

If the absorbance of working reagent is higher than 0.5 at 334 nm the reagent can not be used.

Assay conditions

Wavelength:	340 (334-365) nm
Temperature:	37 °C
Cuvette:	1 cm light path
Method:	endpoint (increasing)
Read against:	reagent blank

- One-reagent procedure

	Blank	Standard	Sample
Dist. water	10 µl		
Standard		10 µl	
Sample			10 µl
Working reagent	1 ml	1 ml	1 ml

Mix and read the absorbance (A) after a 5-minute incubation.

- Two-reagent procedure

	Blank	Standard	Sample
Dist. water	10 µl		
Standard		10 µl	
Sample			10 µl
R1	800 µl	800 µl	800 µl

Mix and wait 1 minute and add:

	Blank	Standard	Sample
R2	200 µl	200 µl	200 µl

Mix and read the absorbance (A) after a 5-minute incubation.

Calibration (37°C, HK test)

S1: Distilled water

S2: Diagnosticum DunaCal Cat. No.: Dcal or

Glucose standard Cat. No.: 450411 or

Roche C.F.A.S. (Calibrator for automated system)

Randox Calibration Serum Level I

Calibration frequency

Two-point calibration is recommended:

- after reagent lot change,

- as required following quality control procedures.

Calculation

$$\frac{A_{\text{sample}}}{A_{\text{standard}}} \times C_{\text{standard}} = C_{\text{sample}}$$

A = Absorbance, C = Concentration

Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P). DunaCont Urine two levels control (DCONU2) is recommended for urine analysis.

PERFORMANCES DATA

The following data were obtained using the Olympus 600 analyzer (37°C). Conversion factor [mmol/l] × 17,97 = [mg/dl]

Linearity

The test is linear between 0,65-33,33 mmol/l (11,7-600 mg/dl).

The limit of detection is 0,018 mmol/l (0,3 mg/dl)

Sensitivity

It is recommended that each laboratory establishes its own range of sensitivity as this is limited by the sensitivity of the spectrophotometer used. Under manual conditions however, a change of 0.001 Abs is equivalent to 0.016 mmol/l (0,29mg/dl) Glucose concentration at 334 nm.

Precision

Reproducibility

n=20	serum		urine			
	sample mmol/l	SD	CV%	sample mmol/l	SD	CV%
	5,8	0,090	1,56	1,24	0,045	3,63
	14,2	0,214	1,51	14,01	0,528	3,77

Repeatability

n=20	serum		
	sample mmol/l	SD	CV%
	4,4	0,065	1,49
	18,1	0,214	1,18

Correlation

Comparative studies were done to compare our reagent with another commercial Glucose assay, on 54 human samples.

The results from these studies are detailed below.

Correlation coefficient: $r=0.9949$

Linear regression: $y \text{ (mmol/l)} = 0.982x + 0.271$

(x=other reagent, y= glucose HK reagent).

Specificity

Bilirubin 855 µmol/l (50 mg/dl), lipid 1000 mg/dl and ascorbic acid 2.84 mmol/l (50 mg/dl) don't interfere with the assay at levels up to the given levels.

Note

Do not use reagents after the expiry date stated on each reagent container label. Do not use products, test solutions and reagents described above for any purpose other than described herein.

For in vitro diagnostic use only.

The following symbols can be used on the labels

In vitro diagnostic device

Batch code

Manufacturer

Catalogue number

CE-marking

This way up

Temperature limitations

Biological risk

Use by (year/month)

Bibliography

Peterson, J.L., Young, D.S., *Anal Biochem.*, 23, (1968); 301
Bondar, R.J.L., Mead, D.C. *Clin. Chem.*, 20, (1974); 586
Young, D.S., Pestaner, L.C., Gibberman, V., *Clin. Chem.*, 5, 10.
Szabó A.: *Klinikai laboratóriumi vizsgálatok és paraméterek (2010)*

Version: 73-DL-2016-14

Date of revision: 2016-12