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

A direct immunoinhibition method for the quantitative determination of high density lipoprotein cholesterol (HDL-C) in serum.

### Principle

Anti human  $\beta$ -lipoprotein antibody in Reagent 1 binds to all lipoproteins of the serum (LDL, VLDL, and chylomicrons) other than HDL. The antigen-antibody complexes formed block enzyme reactions started when Reagent 2 is added. Cholesterol esterase (CHE) and cholesterol oxidase (CO) in Reagent 2 react only with HDL-C. Hydrogen peroxide produced by the enzyme reactions with HDL-C yields a blue color complex upon oxidative condensation of N-ethyl-N-(2-hydroxy-3-sulfo-propyl)-3,5-dimethoxy-4-fluoroaniline, sodium salt (F-DAOS) and 4-aminoantipyrine (4AA) in the presence of peroxidase (POD). By measuring absorbance of the blue color complex produced, at the near optimum wavelength of 593 nm, the HDL-C concentration in the sample can be calculated when compared with the absorbance of the HDL-C Calibrator.

### Reference values

Adult male: >35 mg/dl (>0,9 mmol/l)  
Adult female: >42 mg/dl (>1,1 mmol/l)

It is recommended that each laboratory should assign its own normal range.

### Reagents

1. Pretreatment solution: Store at 2-10°C. Do not freeze.  
30 mmol/l Goods buffer, (pH=7.0), containing 4-AA (0.9 mmol/l), POD (2400 U/l), ascorbate oxidase (2700 U/l), and anti human  $\beta$ -lipoprotein antibody.

2. Enzyme reagent: Store at 2-10°C. Do not freeze.  
30 mmol/l Goods buffer, (pH=7.0), containing CHE (4000 U/l), CO (20000 U/l), and F-DAOS (0.8 mmol/l).

### Safety instructions:

Warning! May Cause allergic skin reaction.

### Samples

Use serum as a specimen. It is recommended to measure HDL-C immediately after collection. Ascorbic acid, bilirubin, and hemoglobin do not have a significant effect on the measurement.

### PROCEDURE

#### Working reagent

Reagents 1 and 2: are ready for use.

#### Stability

without opening: till the expiry date indicated on the label  
after opening: 30 days  
calibration frequency: 30 days  
onboard stability: 30 days  
Stability data are valid only when using new system bottle!

#### Assay conditions

Main wavelength: 600 nm  
Sub wavelength: 700 nm  
Light path: 1 cm  
Temperature: 37 °C  
Method: endpoint (increasing)

#### Pipette into cuvette

	Reagent blank	Standard	Sample
Reagent R1	270 $\mu$ l	270 $\mu$ l	270 $\mu$ l
Standard		3 $\mu$ l	
Sample			3 $\mu$ l
Distilled water	3 $\mu$ l		

Mix and incubate for 5 minutes.

	Reagent R2	90 $\mu$ l	90 $\mu$ l	90 $\mu$ l
Reagent R2	90 $\mu$ l	90 $\mu$ l	90 $\mu$ l	90 $\mu$ l

Mix and incubate for 5 minutes, then read the final absorbance value against reagent blank.

#### Calibration

S1: Distilled water  
S2: Sentinel HC calibrator  
Wako HC calibrator  
Randox HC calibrator

#### Calibration frequency

Two point calibration is recommended:  
- after reagent lot change,  
- as required following quality control procedures.

#### Calculation

$$\frac{A_{\text{Sample}}}{A_{\text{Standard}}} \times C_{\text{Standard}} = C_{\text{Sample}}$$

A = absorbance  
C = concentration

#### Quality control

A quality control program is recommended for all clinical laboratories. The analysis of control material in both the normal and abnormal ranges with each assay is recommended for monitoring the performance of the procedure. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

#### PERFORMANCES DATA

Conversion factor: [mmol/l]=[mg/dl]/38,64

#### Linearity

The test is linear up to 7,15 mmol/l (275 mg/dl).

The reagent is linear between the concentration of 1-180 mg/dl (0.026-4.69 mmol/l). When the concentration of triglyceride in a sample exceeds 1200mg/dl (13.4 mmol/l), dilute the sample with a saline solution, repeat assay and multiply result by the dilution factor.

#### Sensitivity

Absorbance of sample blank is 0.1 or less. Absorbance of a 50 mg/dl (1,30 mmol/l) HDL-cholesterol sample is between 0.07 and 0.34.

#### Specificity

Obtained values of control serum samples with known amount of HDL-cholesterol fall within +/- 10%.

#### Detection Limit

The detection limit for this test is 0,006 mmol/l (0,22 mg/dl)

#### Precision

n=20	Repeatability		
	Average concentration, mg/dl	SD	CV %
1. serum	37,93	0,22	0,58
2. serum	69,79	0,48	0,69

n=20	Reproducibility		
	Average concentration, mg/dl	SD	CV %
1. serum	34,59	0,32	0,93
2. serum	71,3	0,41	0,58
3. serum	92,45	0,59	0,64

#### Correlation

A comparative study has been performed between our reagent and other commercial HDL-cholesterol reagent on 50 human serum samples.

The parameters of linear regression are as follows:

	Serum
n=	50
Regression analysis	y= 0.96x+ 0.065 mmol/l
Correlation coefficient	r=0.998

#### Note

Do not use the pretreatment reagent, which was frozen by mistake.

#### For in vitro diagnostic use only.

The following symbols can be used on the labels

	In vitro diagnostic device		Batch code
	Manufacturer		Catalogue number
	CE-marking		This way up
	Temperature limitations		Biological risk
	Use by (year/month)		

#### Bibliography

- Rifai, N., Warnick, GR. (eds.) *Laboratory Measurement of Lipids, Lipoproteins and Apolipoproteins*. AACC Press, Washington, D.C., USA, 1994.
- Burtis, CA., Ashwood, E.R. (eds.) *Tietz Textbook of Clinical Chemistry*. 2<sup>nd</sup> ed. Saunders, Philadelphia, 1994.
- Gordon, T., Castelli, W.P., Hjortland, M.C., et al., *Am. J. Med.* 62,707-714 (1977)
- Szabó A.: *Klinikai laboratóriumi vizsgálatok és paraméterek (2010)* (ISBN 978-963-9879-75-1)

Version: 76-DL-2017-07

Date of revision: 2017-01