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
Reagent kit for quantitative in vitro determination of gamma-glutamyltransferase (γ -GT) in serum and plasma.

Summary

γ -GT plays an important role in amino acid transport in the course of glutathione metabolism. The enzyme present in the serum is mainly of hepato-biliary origin. Increased enzyme activities are found in association with chronic alcoholism, different toxic liver damages, intra- and extrahepatic cholestasis, acute viral hepatitis, pancreatitis, neoplastic diseases of the liver and pancreas, myocardial infarction as well as with diabetes mellitus.

Method

Optimized IFCC/Szasz method

Reagents

Composition and concentrations

Reagent 1	
Tris buffer, pH: 8.25	125 mmol/l
Glycylglycine	420 mmol/l
Reagent 2	
L- γ -glutamyl-3-carboxy-4-nitroanilide	66 mmol/l

Storage and stability

The reagent is stable up to the end of the indicated month of expiry without opening, if stored at 2 – 8°C, protected from light and contamination is avoided. Do not freeze! Onboard stability after opening and the frequency of calibration is 72 days. The absorbance at 410 nm should not be higher than 1,5

Warnings and precautions

Do not use reagents after the expiry date stated on each reagent container label.

Chemical safety

This product is not classified as dangerous. Safety data sheet is available upon request. The product contains sodium azide. Sodium azide can react with copper and lead plumbing to form explosive metal azides. If disposal into a drain is in compliance with federal, state, and local requirements, flush reagents with a large amount of water to prevent the buildup of azides.

Preparation

The reagent is ready for use.

Sample

Serum, EDTA, heparin or citrate plasma.
Stability in serum
1 month at 2-8°C
5 days at 15-25°C
1 year at -20°C

Expected values and reference range

Serum
Male: 12-52 U/l (0,2 – 0,87 μ kat/l)
Female: 8-34 U/l (0,13 – 0,57 μ kat/l)

It is recommended that each laboratory should assign its own normal range.

Assay procedure

This kit was developed for using in automated systems. Applications can be downloaded from www.diagnosticum.hu or contact us: applikacio@diagnosticum.hu

Wavelength:	410 nm / 700 nm (primary/secondary)
Optical path:	1 cm
Temperature:	37°C
Measurement:	against water blank
Reaction:	kinetic, increasing

	blank	sample or standard
reagent 1	240 μ l	240 μ l
dist. water (diluent)	480 μ l	480 μ l
dist. water (blank)	90 μ l	-
sample or standard	-	90 μ l
	Mix and incubate for 1 minute	
reagent 2	60 μ l	60 μ l
dist. water (diluent)	120 μ l	120 μ l
	Mix and incubate for 1 minute then continuously read the absorbances for 3 minutes	

Calculation

$$\gamma\text{-GT}[\text{U/l}] = \Delta A \text{ sample} / \Delta A \text{ standard} \times \text{standard concentration}[\text{U/l}]$$

Conversion factor

$$[\text{U/l}] = [\mu\text{kat/l}] \times 60$$

Calibration and quality control

- S1: Distilled water
- S2: Diagnosticum DunaCal Cat. No.: Dcal or Roche C.F.A.S. (Calibrator for automated system) or Randox Calibration Serum Level I or

Randox Calibration Serum Level II

Calibration is recommended:

- after opening new reagent batch
- after system maintenance or troubleshooting

For internal quality control, two levels controls are recommended (normal and pathological) at least once a day. The measured values must in the range which was given by the control's manufacturer. Each laboratory should establish corrective measures to be taken if values fall outside the limits. Recommended controls: Diagnosticum DunaCont N (Cat. No.: Dcon-N) and DunaCont P (Cat. No.: Dcon-P)

Performance characteristics

Measuring range

The method is linear in the range 7 – 700 U/l (0,12 – 11,67 μ kat/l)

Interferences

No significant interference was observed by bilirubin up to 1000 μ mol/l bilirubin, triglycerides up to 12 mmol/l, hemoglobin up to 1 g/l and ascorbate up to 4 g/l. Significant interference: >10%.

Limit of detection

The limit of detection is 0,62 U/l

Precision

Repeatability	mean	SD	CV
n = 20	[U/l]	[U/l]	[%]
normal sample	53,1	0,91	1,72
pathological sample	171	2,1	1,23
Reproduceability	mean	SD	CV
n = 10	[U/l]	[U/l]	[%]
normal sample	52,9	1,73	3,27
pathological sample	170,9	2,95	1,73

Method comparison

Comparison with the non-concentrated reagent.

analyser: Advia 2400

number of samples: 150

range: 9 – 1681 U/l

correlation coefficient: 0,9997

regression line equation: $y = 1,03x + 0,789$

(x= normal reagent, y= concentrated reagent)

For in vitro diagnostic use only!

The following symbols can be used on the labels

-  In vitro diagnostic device
-  Manufacturer
-  CE-marking
-  Temperature limitations
-  Use by (year/month)
-  Batch code
-  Catalogue number
-  This way up
-  Biological risk

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