

Rendelési információk:

Katalógusszám Kiszereles
520715200 1 x 5 ml

Felhasználás:

Összefehérje (vizelet) reagens kalibrálására szolgáló stabil folyékony vizeletfehérje standard.

Felhasználásra kész.

A standard koncentrációja LOT specifikus.

Visszavezethetőség:

SRM 927c

Stabilitás:

Felbontás nélkül 2-8°C-on tárolva a lejáratí időig.


Felbontás után 2-8°C-on tárolva minimum 4 hétig.

Figyelmeztetések:

Csak in vitro diagnosztikai célokra használható.


Ne használja a standardet a lejáratí dátum eltelte után!


A címkéken a következő szimbólumok lehetnek

 In vitro diagnosztikai orvostechnikai eszköz

 Gyártó


 CE-jelölés

 Tárolási hőmérséklet

 Lejáratí idő (év/hónap)

 Gyártási szám

 Kódszám

 Tárolási irány

 Biológiai veszély

Ordering information:

Catalogue number Size
520715200,520715N 1 x 5 ml

Preparation:

Total protein ultrasensitive standard for the calibration of total protein ultrasensitive reagent.

Ready to use.

Concentration is LOT specific.

Traceability:

SRM 927c

Stability:

Without opening at 2-8°C until the expiry date.


After first opening at least 4 weeks at 2-8°C.


Warnings:


For in vitro diagnostic use only.


Do not use after the expiry date!


The following symbols can be used on the labels


 In vitro diagnostic device

 Manufacturer

 CE-marking


 Temperature limitations

 Use by (year/month)

 Batch code

 Catalogue number

 This way up

 Biological risk

Verzió/Version: 32B-2015-06

Felülvizsgálat dátuma/Date of revision : 2016-03